



GYNAECOLOGY

 Address: Loujin Industrial park, Shuofang, Wuxi, 214143, P.R.China

 Tel: 86 510 8273 3513 Fax: 86 510 8273 0717


 E-mail: info@gcmedica.com

 Website: www.gcmedica.com

UTERINE ASPIRATION TUBE

FEATURES:

- The tubing is thick enough to eliminate collapse under high vacuum pressure
- 360° rotating handle allows a complete evacuation of the uterine cavity and prevents tubing twisting
- A thumb-indexed ring slides easily over the suction vents for comfortable and positive control during the procedure



Product Code	Description	Qty
GC7050A01	Handle only	50/carton
GC7051C02	2m Tubing with handle, adaptor and flexible cuff	20/carton
GC7051C03	3m Tubing with handle, adaptor and flexible cuff	20/carton
GC7051D02	2m Tubing	20/carton
GC7051D03	3m Tubing	20/carton

Accessories



Product Code	Description	Qty
GC7050A01	Handle for 14mm & 16mm curette, reusable	50/carton

RIGID CURETTE

FEATURES:

- Made of clear rigid plastic, allows for visualization of insertion depth
- Moulded markings for size identification
- Curved and straight design for options
- Smooth tip



Product Code	Description	Qty
GC7049A07	Curved, 7mm	50/carton
GC7049A08	Curved, 8mm	50/carton
GC7049A09	Curved, 9mm	50/carton
GC7049A10	Curved, 10mm	50/carton
GC7049A11	Curved, 11mm	50/carton
GC7049A12	Curved, 12mm	50/carton
GC7049E14	Curved, 14mm	50/carton
GC7049E16	Curved, 16mm	50/carton



Product Code	Description	Qty
GC7049B07	Straight, 7mm	50/carton
GC7049B08	Straight, 8mm	50/carton
GC7049B09	Straight, 9mm	50/carton
GC7049B10	Straight, 10mm	50/carton
GC7049B11	Straight, 11mm	50/carton
GC7049B12	Straight, 12mm	50/carton
GC7049D14	Straight, 14mm	50/carton
GC7049D16	Straight, 16mm	50/carton

FLEXIBLE CURETTE

FEATURES:

- Flexible curette reduce perforation risk and improve patient comfort



Product Code	Description	Qty
GC7049C04	Flexible, 4mm	50/carton
GC7049C05	Flexible, 5mm	50/carton
GC7049C06	Flexible, 6mm	50/carton
GC7049C07	Flexible, 7mm	50/carton
GC7049C08	Flexible, 8mm	50/carton
GC7049C09	Flexible, 9mm	50/carton
GC7049C10	Flexible, 10mm	50/carton
GC7049C11	Flexible, 11mm	50/carton
GC7049C12	Flexible, 12mm	50/carton

UTERINE ASPIRATION KIT



Product Code	Description	Qty
GC7051S02A07	2m Tubing with handle, adaptor, flexible cuff and straight, 7mm rigid curette	10/carton
GC7051S02A08	2m Tubing with handle, adaptor, flexible cuff and straight, 8mm rigid curette	10/carton
GC7051S02A09	2m Tubing with handle, adaptor, flexible cuff and straight, 9mm rigid curette	10/carton
GC7051S02A10	2m Tubing with handle, adaptor, flexible cuff and straight, 10mm rigid curette	10/carton
GC7051S02A11	2m Tubing with handle, adaptor, flexible cuff and straight, 11mm rigid curette	10/carton
GC7051S02A12	2m Tubing with handle, adaptor, flexible cuff and straight, 12mm rigid curette	10/carton
GC7051S02D14	2m Tubing with handle, adaptor, flexible cuff and straight, 14mm rigid curette	10/carton
GC7051S02D14	2m Tubing with handle, adaptor, flexible cuff and straight, 16mm rigid curette	10/carton




Product Code	Description	Qty
GC7051S02B07	2m Tubing with handle, adaptor, flexible cuff and curved, 7mm rigid curette	10/carton
GC7051S02B08	2m Tubing with handle, adaptor, flexible cuff and curved, 8mm rigid curette	10/carton
GC7051S02B09	2m Tubing with handle, adaptor, flexible cuff and curved, 9mm rigid curette	10/carton
GC7051S02B10	2m Tubing with handle, adaptor, flexible cuff and curved, 10mm rigid curette	10/carton
GC7051S02B11	2m Tubing with handle, adaptor, flexible cuff and curved, 11mm rigid curette	10/carton
GC7051S02B12	2m Tubing with handle, adaptor, flexible cuff and curved, 12mm rigid curette	10/carton
GC7051S02E14	2m Tubing with handle, adaptor, flexible cuff and curved, 14mm rigid curette	10/carton
GC7051S02E16	2m Tubing with handle, adaptor, flexible cuff and curved, 16mm rigid curette	10/carton

HSG CATHETER

FEATURES:

- Intended for use in the injection of contrast material in the examination of the uterus and fallopian tubes
- The balloon forms a seal at the internal cervical os to hold the saline or contrast in the uterus
- Rigid insertion sheath
- Radiopaque material of catheter allow clear view during procedure
- Include syringe
- Latex Free



Product Code	Description	Qty
GC7005A11	FR5	10/box
GC7005A12	FR7	10/box

AMNIO HOOK

FEATURES:

- Used during vaginal deliveries to puncture the amniotic sac
- Disposable to prevent infection



Product Code	Qty
GC7053A01	100/carton

IUD REMOVAL HOOK

FEATURES:

- For easy and effective retrieval and removal of an IUD
- Disposable to prevent infection



Product Code	Qty
GC7057A01	100/carton

VAGINAL SPECULUM

FEATURES:

- Crystal clear blade for accurate gynecologic checking
- Unbreakable design available for patient's safety
- Smooth and rounded end for patient's safety and comfort
- With middle screw design
- Color coded



Product Code	Description	Qty	Color
GC7048S01	Size S	40/carton	Red
GC7048S02	Size S, unbreakable	40/carton	Red
GC7048M01	Size M	40/carton	Blue
GC7048M02	Size M, unbreakable	40/carton	Blue
GC7048L01	Size L	40/carton	Green
GC7048L02	Size L, unbreakable	40/carton	Green

FEATURES:

- French type
- Crystal clear blade for accurate gynecologic checking
- With middle screw design



Product Code	Description	Qty	Color
GC7047XS01	Size XS, French type	100/carton	White
GC7047S01	Size S, French type	100/carton	Blue
GC7047M01	Size M, French type	100/carton	Red
GC7047L01	Size L, French type	100/carton	Green

VAGINAL SPECULUM WITH SMOKE EXTRACTOR

FEATURES:

- Sterile Single Use Vaginal Speculum with smoke evacuation tube.
- Smooth round edges increase patient comfort.
- Colour efficiently reflects light, ensuring maximum visibility for the procedure.
- Improved configuration of holes which removes smoke effectively.
- 100pcs per carton.



Product Code	Description	Size
GC7049M01	Vaginal Speculum with Smoke Extractor	M
GC7049L01	Vaginal Speculum with Smoke Extractor	L
GC7049XL01	Vaginal Speculum with Smoke Extractor	XL

